



ARBOR
ASSAYS

DetectX[®]
Superoxide Dismutase
(SOD)
Colorimetric Activity Kit

2 Plate Kit

Catalog Number K028-H1

SPECIES INDEPENDENT

Sample Types Validated:

**Serum, Plasma, Cells, Tissues
and Erythrocyte Lysates**

**Please read this insert completely prior to using the
product.**

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

www.ArborAssays.com


WEB INSERT
TABLE OF CONTENTS

2011-07-01

Background	3
Assay Principle	4
Related Products	4
Supplied Components	4
Storage Instructions	4
Other Materials Required	5
Precautions	5
Sample Types	5
Sample Preparation	5
Standard Preparation	6
Assay Protocol	6
Calculation of Results	7
Typical Data	7-8
Validation Data Sensitivity, Linearity, etc.	8-10
Sample Values	10
Warranty & Contact Information	11
Plate Layout Sheet	12



Short-lived and highly reactive oxygen species (ROS) such as O_2^- (superoxide), $\cdot OH$ (hydroxyl radical), and H_2O_2 (hydrogen peroxide) are continuously generated *in vivo*. In the resting state, the balance between antioxidants and oxidants is sufficient to prevent the disruption of normal physiologic functions; however, either increases in oxidants or decreases in antioxidants can disrupt this balance giving rise to elevated levels of reactive oxygen species (ROS)^{1,2}.



The cellular levels of ROS are controlled by antioxidant enzymes and small molecule antioxidants. The major antioxidant enzymes, superoxide dismutases (SODs), including copper-zinc superoxide dismutase (Cu/ZnSOD), manganese superoxide dismutase (MnSOD) and extracellular superoxide dismutase (EC-SOD), all play critical roles in scavenging O_2^- . Decreased SOD activity results in elevated level of superoxide which in turn leads to decreased NO but increased peroxynitrite concentrations. The major intracellular SOD is a 32-kD copper and zinc containing homodimer (Cu/Zn SOD). The mitochondrial SOD (MnSOD) is a manganese-containing 93-kD homotetramer that is synthesized in the cytoplasm and translocated to the inner matrix of mitochondria. EC-SOD is the primary extracellular SOD enzyme and is highly expressed in many organs. Increased SOD activity levels are seen in Downs Syndrome³ while decreased activity is seen in diabetes, Alzheimer's disease, rheumatoid arthritis, Parkinson's disease, uremic anemia, atherosclerosis, some cancers, and thyroid dysfunction³⁻⁸.

1. Liocher, SI and Fridovich, I., "The Effects of Superoxide Dismutase on H_2O_2 Formation", Free Rad. Biol. Med., 2007, 42:1465-1469.
2. Imlay, JA. "Cellular Defences Against Superoxide and Hydrogen peroxide", Ann Rev. Biochem., 2008, 77:755-776.
3. Torsdottir, G. et al. "Case-control studies on ceruloplasmin and superoxide dismutase (SOD1) in neurodegenerative diseases: a short review". J.Neurol.Sci. 2010. 299(1-2):51-54.
4. Giacco, F & Brownlee, M. "Oxidative stress and diabetic complications". Circ.Res. 2010. 107(9):1058-1070.
5. Bae, S-C, et al. "Inadequate antioxidant nutrient intake and altered plasma antioxidant status of rheumatoid arthritis patients", J.Amer.Coll.Nutr. 2003. 22(4):311-315.
6. Akbostanci, MC, et al. "Erythrocyte superoxide dismutase activity differs in clinical subgroups of Parkinson's disease patients". Acta Neurol.Belg. 2001. 101:180-183.
7. Shainkin-Kestenbaum, R, et al. "Reduced superoxide dismutase activity in erythrocytes of dialysis patients: a possible factor in the etiology of uremic anemia". Nephron. 1990. 55(3):251-253.
8. Saito, T. "Superoxide dismutase level in human erythrocytes and its clinical application to patients with cancers and thyroidal dysfunctions". Hokkaido Igaku Zasshi. 1987. 62(2):257-268.
9. Okado-Matsumoto, A, and Fridovich, I., "Subcellular Distribution of Superoxide Dismutases (SOD) in Rat Liver". J. Biol. Chem. 2001, 276:38,388-38,393.

WEB INSERT ASSAY PRINCIPLE

2011-07-01

The DetectX[®] Superoxide Dismutase (SOD) Activity Kit is designed to quantitatively measure SOD activity in a variety of samples. The assay measures all types of SOD activity, including Cu/Zn, Mn, and FeSOD types. Please read the complete kit insert before performing this assay. A bovine erythrocyte SOD standard is provided to generate a standard curve for the assay and all samples should be read off of the standard curve. Samples are diluted in our specially colored Sample Diluent and added to the wells. The Substrate is added followed by Xanthine Oxidase Reagent and incubated at room temperature for 20 minutes. The Xanthine Oxidase generates superoxide in the presence of oxygen, which converts a colorless substrate in the Detection Reagent into a yellow colored product. The colored product is read at 450 nm. Increasing levels of SOD in the samples causes a decrease in superoxide concentration and a reduction in yellow product. The activity of the SOD in the sample is calculated after making a suitable correction for any dilution, using software available with most plate readers. The results are expressed in terms of units of SOD activity per mL.

RELATED PRODUCTS

DetectX[®] Kits

Hemoglobin Dual Range Detection Kit	Catalog Number K013-H1
Glutathione Fluorescent Detection Kits	Catalog Number K006-F1/F5
Glutathione Colorimetric Detection Kit	Catalog Number K006-H1
Glutathione S-Transferase Fluorescent Activity Kit	Catalog Number K008-F1
Glutathione Reductase Fluorescent Activity Kit	Catalog Number K009-F1
Nitric Oxide Colorimetric Detection Kit	Catalog Number K024-H1



WEB INSERT
SUPPLIED COMPONENTS

2011-07-01

Clear 96 well Half Area Plates	2 Plates	Catalog Number X018-2EA
Superoxide Dismutase Standard	1 vial	Catalog Number C098-1EA
1 Unit/vial of bovine Erythrocyte Superoxide Dismutase lyophilized.		
Assay Buffer	50 mL	Catalog Number X100-50ML
Buffer containing detergents, stabilizers and dye.		
Xanthine Oxidase Buffer	6 mL	Catalog Number X102-6ML
Buffer containing detergents and stabilizers.		
Xanthine Oxidase Concentrate	225 μ L	Catalog Number C099-225UL
A 25X concentrated suspension of Xanthine Oxidase.		
Substrate Diluent	12 mL	Catalog Number X101-12ML
Substrate buffer.		
Substrate Concentrate	1.1 mL	Catalog Number C100-1.1ML
A 10X concentrate of the Detection Reagent.		

STORAGE INSTRUCTIONS

All components of this kit should be stored at 4°C until the expiration date of the kit.

Once reconstituted, the Superoxide Dismutase Standard should be aliquoted and stored at -20°C.

WEB INSERT
OTHER MATERIALS REQUIRED

2011-07-01

2 mM Potassium Cyanide solution for inhibition of Cu/Zn and extracellular SOD.

Repeater pipet with disposable tips capable of dispensing 25 and 50 μ L.

96 well plate reader capable of reading optical absorption at 450 nm.

Software for converting optical density (OD) readings from the plate reader and carrying out four parameter logistic curve (4PLC) fitting. Contact your plate reader manufacturer for details.

PRECAUTIONS

As with all such products, this kit should only be used by qualified personnel who have had laboratory safety instruction. The complete insert should be read and understood before attempting to use the product.

SAMPLE TYPES AND PREPARATION

Samples that need to be stored after collection should be stored at -70°C or lower, preferably after being frozen in liquid nitrogen. This assay has been validated for serum, plasma and erythrocyte lysates. Samples containing visible particulate should be centrifuged prior to using.

Some serum and plasma samples may contain significant hemoglobin concentrations and the optical density at 450nm determined prior to running the assay. After addition of the Substrate solution to all the used wells the optical density at 450 nm should be read and subtracted from the optical density recorded at the end of the 20 minute incubation.

In order to measure cytosolic and/or mitochondrial SOD from samples, the supernatant samples prepared on page 7 should be centrifuged at $10,000 \times g$ for 15 minutes at 4°C . The supernatant will contain the cytosolic SOD and the cell pellet will contain mitochondrial SOD. Process the cell pellet as described for Cell Lysates.

To determine Mn SOD activity treat the sample with 2 mM potassium cyanide. Addition of cyanide will inactivate other SOD enzymes.



Samples should be kept on ice to maintain enzyme activity.

Cell Suspensions and Adherent Cells

1. Centrifuge $> 1 \times 10^6$ cells in suspension at 250 x g for 10 minutes at 4°C. Discard the supernatant.
2. Resuspend the cell pellet in ice-cold PBS and transfer to a microtube on ice. Centrifuge, discard supernatant, and place on ice.
3. Wash $> 1 \times 10^6$ adherent cells with PBS prior to being harvested by gentle trypsinization. Transfer to a tube on ice and centrifuge at 250 x g for 10 minutes at 4°C and discard the supernatant.
4. Wash the pellet with ice-cold PBS and centrifuge at 250 x g for 10 minutes at 4°C.
5. Homogenize or sonicate the pellet in 0.5-1 mL of PBS per 100 mg of cells. Centrifuge at 1,500 x g for 10 minutes at 4°C.
6. Collect the supernatant and assay immediately, or store at -80°C. Dilute in Assay Buffer prior to measuring SOD activity.

Tissue Samples

1. Wash tissue thoroughly with ice cold PBS prior to processing.
2. Homogenize as described in steps 5 and 6 for Cell Suspensions above.

Plasma Samples

1. Collect plasma in heparin or EDTA tubes. Centrifuge the sample at 700 x g for 15 minutes at 4°C. Aspirate off the pale yellow supernatant. Assay immediately or freeze at -80°C.
2. Plasma should be diluted at least 1:5 by taking one part of plasma and adding 4 parts of Assay Buffer prior to assaying. Further dilutions in Assay Buffer may be needed.

RBC/Erythrocytes

1. Erythrocytes can be lysed by taking the pelleted RBCs from the Plasma step above and adding 4 volumes of ice cold deionized water.
2. Centrifuge at 10,000 x g for 15 minutes at 4°C to remove debris.
3. Lysed RBCs will exhibit high background color. After adding the Detection Substrate solution read the blank OD at 450 nm prior to addition of Xanthine Oxidase reagent.
4. RBCs should be diluted at least 1:100 prior to running in the assay. Take 10 μ L of lysed RBCs and add to 990 μ L of Assay Buffer prior to assaying. Further dilutions in Assay Buffer may be needed.

Serum Samples

1. Collect blood in serum tubes. Centrifuge the sample at 700 x g for 15 minutes at 4°C. Aspirate off the serum supernatant.
2. Serum should be diluted at least 1:5 with the provided Assay Buffer prior to assaying. Take one part of serum and adding 4 parts of Assay Buffer prior to assaying. Further dilutions in Assay Buffer may be needed.


WEB INSERT STANDARD PREPARATION

2011-07-01

Standard Preparation

Superoxide Dismutase Standards are prepared by labeling six tubes as #2 through #7. The lyophilized vial of SOD Standard is used for Standard 1. Add 250 μL of Assay Buffer to the vial, vortex and let stand at room temperature for 5 minutes. Pipet 75 μL of Assay Buffer into tubes #2 to #7. Carefully add 75 μL of the SOD Standard from the vial to tube #2 and vortex completely. Repeat this for tubes #2 through #7. The activity of SOD in the SOD Standard vial and tubes 2 through 7 will be 4, 2, 1, 0.5, 0.25, 0.125 and 0.0625 U/mL.

Use all Standards within 2 hours of preparation. Aliquot the reconstituted SOD vial and freeze at -20°C .



	Std 1	Std 2	Std 3	Std 4	Std 5	Std 6	Std 7
Assay Buffer Vol (μL)	250	75	75	75	75	75	75
Addition	Stock	Std 1	Std 2	Std 3	Std 4	Std 5	Std 6
Vol of Addition (μL)	-	75	75	75	75	75	75
Final Activity (U/mL)	4	2	1	0.5	0.25	0.125	0.0625

Xanthine Oxidase Preparation

Vortex the suspension of Xanthine Oxidase prior to pipetting. Pipet from the base of the tube.

	1/2 Plate	1 Plate	1.5 Plates	2 Plates
Xanthine Oxidase	50 μL	100 μL	150 μL	200 μL
Xanthine Oxidase Buffer	1.2 mL	2.4 mL	3.6 mL	4.8 mL
Final Mixture	1.25 mL	2.5 mL	3.75 mL	5 mL

Substrate Preparation

Vortex the vial of Concentrate prior to pipetting.

	1/2 Plate	1 Plates	1.5 Plates	2 Plates
Substrate Concentrate	250 μL	500 μL	750 μL	1 mL
Substrate Diluent	2.25 mL	4.5 mL	6.75 mL	9 mL
Final Mixture	2.5 mL	5 mL	7.5 mL	10 mL

WEB INSERT ASSAY PROTOCOL

2011-07-01

Use the plate layout sheet on the back page to aid in proper sample and standard identification.

1. Pipet 10 μ L of samples or appropriate standards into duplicate wells in the plate.
3. Pipet 10 μ L of Assay Buffer into duplicate wells as the Zero standard.
4. Add 50 μ L of the Substrate Preparation to each well using a repeater pipet.

NOTE: If your samples have significant yellow coloration then pre-read the optical density at 450 nm.

5. Add 25 μ L of the Xanthine Oxidase Preparation to each well using a repeater pipet.
6. Incubate at room temperature for 20 minutes.
7. Read the optical density at 450 nm.

WEB INSERT
CALCULATION OF RESULTS

2011-07-01

Average the duplicate OD readings for each standard and sample.

If your sample was visibly colored, then subtract the pre-read optical density at 450 nm from the subsequent Xanthine Oxidase reaction optical density after 20 minutes.

Inhibition values are sometimes quoted for SOD activity. Inhibition values can be obtained by dividing the measured Mean OD for the standard or samples by the Mean OD for the Zero standard (No SOD) and multiplying the result by 100. Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic curve (4PLC) fit. The sample activities obtained should be multiplied by the dilution factor to obtain neat sample values.


TYPICAL DATA

Sample	Mean OD	% Inhibition	SOD Activity (U/mL)
Standard 1	0.055	8.8	4
Standard 2	0.085	13.6	2
Standard 3	0.157	25.0	1
Standard 4	0.235	37.5	0.5
Standard 5	0.362	57.7	0.25
Standard 6	0.479	76.4	0.125
Standard 7	0.563	89.8	0.0625
Zero	0.627	100	0
Sample 1	0.297	47.3	0.360
Sample 2	0.122	19.4	1.307

**Always run your own standard curves for calculation of results.
Do not use these data.**



Typical Standard Curve



Always run your own standard curves for calculation of results.
Do not use these data.

VALIDATION DATA

Sensitivity and Limit of Detection


Sensitivity was calculated by comparing the ODs for twenty wells run for each of the zero and standard #7. The detection limit was determined at two (2) standard deviations from the zero along the standard curve.

Sensitivity was determined as 0.044 U/mL.

Linearity

Linearity was determined by taking two samples, one with a high known SOD activity and the other with a lower SOD activity and mixing them in the ratios given below. The measured activities were compared to the expected values based on the ratios used.

High Sample	Low Sample	Observed Activity (U/mL)	Expected Activity (U/mL)	% Recovery
80%	20%	0.920	0.987	93.2
60%	40%	0.776	0.789	98.4
40%	60%	0.534	0.590	90.4
20%	80%	0.414	0.392	105.6
			Mean Recovery	96.9%



Intra Assay Precision

Three samples diluted in Assay Buffer were run in replicates of 20 in an assay. The mean and precision of the calculated concentrations were:

Sample	SOD Activity (U/mL)	%CV
1	0.407	4.6
2	0.726	7.3
3	1.203	16.8

Inter Assay Precision

Three samples diluted in Assay Buffer were run in duplicates in sixteen assays run over multiple days by four operators. The mean and precision of the calculated concentrations were:

Sample	SOD Activity (U/mL)	%CV
1	0.356	10.5
2	0.653	6.1
3	1.277	13.8

WEB INSERT
SAMPLE VALUES

2011-07-01

Five random adult human serum and plasma samples were diluted in Assay Buffer and run in the assay. The serum samples ranged from 1.95 to 4.60 U/mL with an average of 3.44 U/mL. EDTA plasma samples ranged from 2.24 to 3.56 U/mL with an average 2.89 U/mL. Five samples of RBC from EDTA were normalized to hemoglobin levels using the DetectX[®] Hemoglobin Detection kit, K013-H1. The RBC activities ranged from 748 to 1,507 U/g Hgb with an average of 1052 U/g Hgb.

WEB INSERT
LIMITED WARRANTY

2011-07-01

Arbor Assays warrants that at the time of shipment this product is free from defects in materials and workmanship. This warranty is in lieu of any other warranty expressed or implied, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose.

We must be notified of any breach of this warranty within 48 hours of receipt of the product. No claim shall be honored if we are not notified within this time period, or if the product has been stored in any way other than outlined in this publication. The sole and exclusive remedy of the customer for any liability based upon this warranty is limited to the replacement of the product, or refund of the invoice price of the goods.

CONTACT INFORMATION

For details concerning this kit or to order any of our products please contact us:

Arbor Assays

1514 Eisenhower Place
Ann Arbor, Michigan 48108 USA
Phone: 734-677-1774
Fax: 734-677-6860
Web: www.ArborAssays.com

E Mail Addresses:

Info@ArborAssays.com
Orders@ArborAssays.com
Technical@ArborAssays.com
Contracts@ArborAssays.com

DetectX®, ThioStar® and the Arbor Assays logo are all registered trademarks.



	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

