

# VEQ 360°



**Una risorsa indispensabile per confermare e migliorare la qualità di tutto il tuo laboratorio.**

Catalogo 2024  **INSTAND**  
programmi dedicati al Complemento

Estratto dal catalogo originale instand 2024

## Complement Analysis 01 – Auto-anti-C1q

Prohaszka / Kirschfink

245

Samples: 2 samples of 300 µl Sample properties: Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 12 months		additional sample set not available
Auto-anti C1q (IgG)	Auto-anti C1q (IgG) (qualitative)	
	Mar. (2)	Nov. (6)
Registration deadline:	29.12.23	09.08.24
Sample Shipment:	27.02.24	08.10.24
Closing Date:	29.03.24	08.11.24

## Complement Analysis 02 – Functions and proteins

Prohaszka / Kirschfink

246

Samples: 2 samples of 1 ml Sample properties: Lyophilised serum Shipping Information: Shipping at ambient temperature Certificate Validity: 12 months			Maximum order: 3 sets of samples
Alternative pathway C1q Factor H Lectin pathway	C1 inhibitor: activity C3 Factor I	C1 inhibitor: protein C4 Classical pathway	
Analysis of the total function of the 3 complement activation pathways (classical, alternative and lectin pathways) either by functional ELISA (WIELISA), by haemolytic (CH50, AH50) or by liposome-based assay (classical pathway). Participation with one parameter (e.g. CH50) is possible.			
	Mar. (2)	Nov. (6)	
Registration deadline:	29.12.23	09.08.24	
Sample Shipment:	27.02.24	08.10.24	
Closing Date:	29.03.24	08.11.24	

## Complement Analysis 03 – Complement activation products

Prohaszka / Kirschfink

247

Samples: 2 samples of 300 µl Sample properties: Lyophilised plasma Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months		Maximum order: 3 sets of samples
C3a, C3d, Bb, sC5b-9		
Participation with one parameter (e.g. sC5b-9) possible		
		<b>Nov. (6)</b>
Registration deadline:		09.08.24
Sample Shipment:		08.10.24
Closing Date:		08.11.24

Samples: 2 samples of 300 µl Sample properties: Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months	additional sample set not available
C3-nephritic factor (qualitative)	
<div><div>Registration deadline:</div><div>Sample Shipment:</div><div>Closing Date:</div></div> <div><div>Nov. (6)</div><div>09.08.24</div><div>08.10.24</div><div>08.11.24</div></div>	

Samples: 2 samples of 300 µl Sample properties: Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months		additional sample set not available
Auto-anti FH-(qualitative)	Auto-anti-FH (quantitative)	
<div>Nov. (6)</div> <div>Registration deadline: 09.08.24 Sample Shipment: 08.10.24 Closing Date: 08.11.24</div>		

Samples: 2 samples of 300 µl Sample properties: Serum Shipping Information: Shipping at ambient temperature Certificate Validity: 24 months		additional sample set not available
Auto-anti-C1 inhibitor IgA (qualitative) Auto-anti-C1 inhibitor IgM (qualitative)	Auto-anti-C1 inhibitor IgG (qualitative)	
Nov. (6)		
Registration deadline:	09.08.24	
Sample Shipment:	08.10.24	
Closing Date:	08.11.24	