



## Product Information Pipette

### **Product information for BP10/BP20/BP200/BP1000:**

#### **Introduction**

The Bio Basic Inc. (BBI) pipette is a volumetric instrument designed to measure and transfer liquids precisely and safely. It can measure and transfer volumes from 0.5 - 1000µl depending on the model.

The BBI pipette has a digital counter which displays the volume. The volume is adjusted by turning the orange knurled wheel. The maximum volume of the pipette is shown on the pushbutton and corresponds to the model number.

#### **All BBI pipettors come with a 1-year manufacturer's warranty"**

The manufacturer warrants that the Electronic Pipettor meets its current specifications for 1 year after shipment, provided it is used and cared for in accordance with its instructions. If for any reason the product fails to meet any express warranties, as your sole remedy, Bio Basic will repair or replace the unit or refund the purchase price.

Item No.	Volume Range	Accuracy (%)
BP10	0.5-10ul	± 2.5
BP20	2-20ul	±1.8
BP200	20-200ul	± 0.8
BP1000	100-1000ul	± 0.8

